

Remove crowns & bridges in 2 minutes risk-free painlessly... and reuse them !

Why take risks?

1. WAMkey is a simple to use solution
2. It will not cause any pain to your patients
3. It's quick (usually less than 2 minutes per crown)
4. In most situations, you can reuse the crown provisionally, Thus, you can save a lot of time. In few cases, you can even reuse them permanently.

This kit has everything you need:

3 sizes of WAMkeys No1, No2, No3

3 CMV (Ceramometalivore) autoclavable burs for effortless, vibration free cutting of metal and porcelain

Easy to follow, complete instruction CDRom

All components in this kit are fully autoclavable (130°C/275°F)