Curing Lights: WL-070

WL-070





| The continue of the continue

Rating: Not Rated Yet

1 / 3

Curing Lights: WL-070

Price

\$1.00

Ask a question about this product

Description

Dentmate Curing Light combination cord/cordless use.>800mw/cm2.

Wireless charging stand

No need to plug in the cable; No pogo pin required.

Built-in intensity indicator

Easy to interpret due to the built-in intensity indicator.

Simple operations

Clear display with simple functions, including the standard mode and ramp mode; easy to use.

Cure fast in 10 secs.

Having the high power 5W LED, this device can cure the composites over 3mm in 10 secs.

Never Shutdown After Continuous Use

Having the smart heat management system, the device will not be overheated after continuous use.

Varieties of optical fiber light guide rod

There are 3, 5, and 8 mm black or transparent optical fiber light guide rod to be applied.

Color Choices

(White, Black, Orange, Blue, Green, Red, Purple)

Specification

Handpiece

Light source: High power dental blue 5W LED (average lifespan: 100,000 hours)

Wavelength: 440-480 nm, peak: 460 nm

Radiant Intensity: Up to 1,200 mW/cm²

Light Guide Rod: Ø8mm optical fiber, Sterilizable in autoclave (max. T. 135? for 20 minutes)

Weight: 120 g (with battery & light guide rod)

Dimensions: Ø26 (max) x L156 mm

Wireless charging stand

Weight: 85 g

2 / 3

Curing Lights: WL-070

Dimensions: W45 x L135 x H40 mm

Battery

Type: Lithium-ion, 3.7V/800mAh

Charging Time: About 3 hours

Use Time: 40 minutes (200 times in 10 Sec. mode)

Power Supply

Type: Mircro USB

Input: 100 ~ 240 V, 50-60 Hz

Output: DC 5V/2A

Reviews

There are yet no reviews for this product.

3 / 3