

**Filratch Manual Pin Inserter for soft Dentin. 0.6mm Red**



## Better by design **FILPIN** restoration retention system

FILPIN has been design engineered to be faster and easier to place even in difficult cases, safer to place by reducing the risk of tooth cracking and crazing, and stronger in use via its unique 'retentive' thread.

### When to use

- To enhance the retention of all types of restorations to dentin.
- Compared to using non-bonded amalgam only, FIL PIN can increase the retention of restorations to dentin by around 3 to 5 times.  
*(L.S. et al., J Dent Res '94; 73:187)*
- To stabilise by cross-pinning MOD cavities.
- To prevent rotation of restorations when a single post is used.



easier

**99.9% pure titanium**  
biocompatible, will not corrode, compatible with all dental materials.



**Self-sharps free time, every time**  
once optimum depth of 2mm reached.

faster

**Self-threading, self-aligning design**  
to make placement faster, easier and safer as reduces the risk of cracking and crazing.



**Long and slim**  
for better visibility during placement, also helps pin to self-align.

safer

**Unique 'retentive' thread**  
maintains retention strength without causing the torsional stresses which can lead to cracking or crazing.



**Easy to bend after insertion**  
without breaking the pin or tooth.

stronger

**Pin hole DENT HEAD**  
1.5mm at the shoulder so pin hole is drilled to the exact depth (2mm).  
Colour coded with pin for ease of selection.



Filratch Manual Pin Inserter for soft Dentin. 0.6mm Red

Rating: Not Rated Yet

Price

\$1.00

[Ask a question about this product](#)

### Reviews

There are yet no reviews for this product.